

# EUROSPIN II TEST SYSTEM


## Product Specifications

Item	Code	Description	
<b>INDIVIDUAL ITEMS</b>			
1	Q50011	<b>TO1 – UNIFORMITY</b> A chamber filled with a slightly acidic solution and doped to achieve T1 of 300ms at 8MHz. Evaluates Uniformity, Ghosting, Image Signal to Noise Ratio and Uniformity of Image Signal to Noise. Dimensions: 200mm x 62mm Weight: 2.0kg filled Weight: 0.6kg empty Solution to fill: 1425ml	
2	Q50021	<b>TO2 – SLICE</b> A chamber like TO1 but also containing an array of geometric shapes. Evaluates Geometric Distortion, and Slice Profile and Width. Dimensions: 200mm x 62mm Weight: 2.0kg filled Weight: 0.9kg empty Solution to fill: 1153ml	
3	Q50031	<b>TO3 - SLICE POSITION</b> A chamber like TO1 but also containing a lattice of fine rods. Evaluates Slice Warp and Slice Position (single and multislice). Dimensions: 200mm x 52mm Weight: 1.7kg filled Weight: 0.7kg empty Solution to fill: 1020ml	
4	Q50041	<b>TO4 – RESOLUTION</b> A chamber like TO1 but also containing objects of various precise sizes and profiles. Evaluates Resolution by Bar Pattern and the Modulation Transfer Function. Dimensions: 200mm x 52mm Weight: 1.7kg filled Weight: 0.9kg empty Solution to fill: 832ml	
5	Q50051	<b>TO5 - CONTRAST</b> , with set of 18 Test Tubes and Gels. A chamber with 12 cavities into which vials of unique tissue mimicking gel can be inserted. Evaluates T1 & T2 Precision, Accuracy, and Contrast. Dimensions: 200mm x 52mm Weight: 1.4kg filled Weight: 0.6kg empty Solution to fill: 740ml	
6	Q50071	<b>Large Carrying Case</b> (for TO1 to TO4) Special design with no magnetic parts. Dimensions 54cm x 52cm x 19.2cm Weight: 3.7kg	
7	Q50081	<b>Small Carrying Case</b> (for TO5 and 18 Gel Tubes) Special design with no magnetic parts. Dimensions 33cm x 29cm x 16cm Weight: 1.3kg	

# EUROSPIN II TEST SYSTEM

## Product Specifications

Item	Code	Description
8	Q50061	<b>USER GUIDE</b> Dimensions: A4, 80 pages Weight: 450g
9	Q50311	<b>Top-Up Solution</b> , 2.5 litres Weight: 4.5kg
10	Q50321	<b>Top-Up Solution</b> , 500ml Weight: 0.92kg
11	Q50201	<b>LOADING ANNULUS</b> Dimensions: 240mm x 310mm Weight: 3.8kg filled Weight: 2.2kg empty Solution to fill: 1890ml
12	Q50211	<b>Annulus Carrying Case</b> Dimensions: 41cm x 34cm x 36cm Weight: 4.3kg
13	Q50091	<b>Set of 18 Test Tubes &amp; Gels</b> Each 24mm diameter x 100mm length Weight: 1.1kg
14	Q50101	<b>Individual Test Tube &amp; Gel</b> Weight of individual Tube: 60g



Prepared November 2011.

## DIAGNOSTIC SONAR LTD

Baird Road, Kirkton Campus, LIVINGSTON EH54 7BX

Telephone: 01506-411877, Fax: 01506-412410

Email: [sales@diagnosticsonar.com](mailto:sales@diagnosticsonar.com)

Web: [www.diagnosticsonar.com](http://www.diagnosticsonar.com)